

CRABP1-cellular retinoic acid binding protein 1 Gene

Reference: PTXBC022069

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CRABP; CRABP-I; CRABPI; RBP5; cellular retinoic acid-binding protein 1; cellular retinoic acid-binding protein I; cellular retinoic acid binding protein 1		

For research use only.

GENE DESCRIPTION	cellular retinoic acid binding protein 1
GENE ID	1381
ORIGIN SPECIES	Human
ACCESSIONS	BC022069
VECTOR	pDONR223
NCBI SYMBOL	CRABP1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcccaactcgcggcacctggaagatgcgcagcagcagagaatttcgacgagctgtcaaggcactgggtgtgaacccatgctgaggaaagtggccgtagcggctgtccaagccgcacgtggagatccgagacggggatcagttctacatcaagacatccaccaggtg
cgcaccactgagatcaacttcaagtcgggagaaggctttgaggagagaccgtggacggacgcaagtgcaggagtttagccacttgggagaatgagaacaagatccactgcacgcaaaccttcttgaaggggacggcccaaacctactggaccagtgcagctggccaacgatgaa
cttatcctgacgtttggcgccgatgacgtggctgcaccagaatttatgccgagagtga