

KIF6-kinesin family member 6 Gene

Reference: PTXBC022074

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	kinesin-like protein KIF6; C6orf102; dJ1043E3.1; dJ137F1.4; dJ188D3.1; kinesin family member 6		

For research use only.

GENE DESCRIPTION	kinesin family member 6
GENE ID	221458
ORIGIN SPECIES	Human
ACCESSIONS	BC022074
VECTOR	pDONR223
NCBI SYMBOL	KIF6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

```
atgagagaggaaatgctattagatccaggaggctttgaaatctcaagaggaccgctgacagcgttaccatcgatgacaacaacagattctgaagccaaggccctggagaaagtataaatgaagcaagaagtaaaattggtcacctgaaggaagaa
atcaccacggcgcatafacgcaagtagccctaggaatctggaaaacatggccgtgctctgatgccagaccagcaggagagaagctgcatcaactggagaaagaaagagaggtataaaacaatgtcactcgcctgaaagccctgaaagtgagatcgagcactgcag
ctgctcatggacaagccaaggtgaagctacagaagaggttgaagctctggggcagaggagcccaaacctgcaggtaaattctccagcagtgaaattcactcgtacacagaagccattctccagacatctgactccagcatgaatggtccaactctcttaacaaaagttctgg
aggctgggaagtccaagatcaaggcactggcagattcgtctgtgatggaatgccaggaatactgcccctgctgcccagtcacacagcagaacagagcagcaccagcaccactggaagacagatcccaagagccagtgctcctcctcaccggagac
agccagacggactcggacatcatgccttcaagccagacagagcattctgagaagcaatggttgggaagcaatgga
```