

LYSMD2-LysM, putative peptidoglycan-binding, domain containing 2 Gene

Reference: PTXBC022075

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	lysM and putative peptidoglycan-binding domain-containing protein 2; LysM, putative peptidoglycan-binding, domain containing 2; LysM domain containing 2		

For research use only.

GENE DESCRIPTION	LysM, putative peptidoglycan-binding, domain containing 2
GENE ID	256586
ORIGIN SPECIES	Human
ACCESSIONS	BC022075
VECTOR	pDONR223
NCBI SYMBOL	LYSMD2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaacagattaaaggccaataaactgtttaccaatgattgtatattctgaagaaactttgaacatcccagttatcatcagagaagcctttgtgttaaggacttaactcattgattctccagaaaatgaaactgctgataacagttttctcaggaaggagccagtggtggccgggga
agacctccctccctccagctcctcaagaatctgatgttcagctgtgcagcctgaggaagtgtcagccagagatttctgcagagacttgacttgcagattaagttatcaacagcagccaagaagctaaaagaagagagcagagatgaagaaagtcctatgcaacttcctctatcaca
gtag