

KLHDC9-kelch domain containing 9 Gene

Reference: PTXBC022077

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	KARCA1; kelch domain-containing protein 9; kelch/ankyrin repeat-containing cyclin A1-interacting protein; kelch domain containing 9		

For research use only.

GENE DESCRIPTION	kelch domain containing 9
GENE ID	126823
ORIGIN SPECIES	Human
ACCESSIONS	BC022077
VECTOR	pDONR223
NCBI SYMBOL	KLHDC9
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaacagcttgcaaggctgtgagcagtgaggcaggggtccagaaggggccccatggactacggcatcactcatgttctgtgtcgggccctttgctgtgctgtttggaggagaaactctgaccagagctagagacaccatctgcaatgatctctacatctatgatactcgacatctctc
ctttgtgttccacttccctgtgagatcgtgggatgaaacgcattggccatcgacctgcctttggaatgatcagctttacctggttgggggttttggaggatggcaggacagccagtcacaggtttgcatcctggactttatctaa