

KHDC1-KH homology domain containing 1 Gene

Reference: PTXBC022080

| | | | |
|----------------|--|---------------------|----------------|
| Product type: | DNA | Storage Condition: | -20°C |
| Size/Quantity: | 2 ug | Delivery Condition: | +4°C, Blue Ice |
| Synonyms | C6orf147; C6orf148; Em:AC019205.8; NDG1; bA257K9.4; KH homology domain-containing protein 1; KH homology domain containing 1 | | |

For research use only.

| | |
|-------------------------|---|
| GENE DESCRIPTION | KH homology domain containing 1 |
| GENE ID | 80759 |
| ORIGIN SPECIES | Human |
| ACCESSIONS | BC022080 |
| VECTOR | pDONR223 |
| NCBI SYMBOL | KHDC1 |
| LENGTH | ProteoGenix bp |
| SEQUENCE PRIMERS | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

ORF (Only for reference, product provided without STOP codon)

atggacatgggaacgagtgctctcagcaagaagccgtggtggacctgcctcaaaacttcatgcaccaatggtgttcacatggaagaggaccaggagctcatctcgggcatggtgacacatacctcgcctgcatgaggtgcacagccacacccttattcagctggagagttggttcacagctacagccagactcgtgtgactgfttgaccacacagggcaaggcagtgctgctgcacatgtttgtgtgtgggagccaggactcctatcatcatgctcagggcctggagatgctggaactgtccgaagccagcctctgaccaacgatgacctgtcacctccattagtggtccaccatacctggagacctgtcttggctcccaggataagtggaactgtgtcttagcgttccctcagcctcccctaccaagtgtggtgtgtcgggttccattgagttcactgtatccgttaa