



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

ENHO-energy homeostasis associated Gene

Reference: PTXBC022101

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

C9orf165; UNQ470; adropin; GAAI470; energy homeostasis-associated protein; energy homeostasis associated

For research use only.

GENE DESCRIPTION	energy homeostasis associated
GENE ID	375704
ORIGIN SPECIES	Human
ACCESSIONS	BC022101
VECTOR	pDONR223
NCBI SYMBOL	ENHO
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggggggcagccalcctccagggggccctcatgccatcgtcgcaacggctctgtggggctcttgcgtctcgtctctgggtcatcctctgtgggcctgcattctcgtctgccgacgttgactctctctgaatccagaactcaagatgcgcctgcttgaggctggaccacggcgccc