

AFTPH-aftiphilin Gene

Reference: PTXBC022247

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	Nbla10388; aftiphilin; aftiphilin protein		
For research use only.			

GENE DESCRIPTION	aftiphilin
GENE ID	54812
ORIGIN SPECIES	Human
ACCESSIONS	BC022247
VECTOR	pDONR223
NCBI SYMBOL	AFTPH
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtatgcagcaggattgggtatgttagagcccaccaaggaaccactgaaaccactttctgtcagaaaaatagcttcacgtcagacagccaccatgtcaccagatatgaacacatgtacatctgatcagttccaggagtctctaccaccgtccagtttgactggagtagcagtgcc
cttactaaccttttagatgggtggtcgaggtgtgatgagtaacaacttctaagctggaaatctccacctcaagcctcaagtgactgatgcatttgcagactcatgtctacagtagagaagacaagcacatctaccaggaaaccgaaaagagaaggccacctaagtgaagaagctatca
agggtgatcgtggccttcctgacttaacattcatgcatgccaagggtgtga