

ZFAND1-zinc finger, AN1-type domain 1 Gene

Reference: PTXBC022251

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AN1-type zinc finger protein 1; zinc finger, AN1-type domain 1; zinc finger AN1-type containing 1		

For research use only.

GENE DESCRIPTION	zinc finger, AN1-type domain 1
GENE ID	79752
ORIGIN SPECIES	Human
ACCESSIONS	BC022251
VECTOR	pDONR223
NCBI SYMBOL	ZFAND1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggagttggacatgggcagcactgccaggtggagcattgccggcagcagatttcttccatttctgtgtgatgtgtcaggaatatttgccttgaacacagaagcagggagtcctggtgtcctgaggtgactgtaataatgagagactgaagacagatcaacatacatctta
cccatgctctttcaagactgtgctgagagagaacttggcagttatgtcttattgtgagaagaatttgcctgagacaccgtcatcagtcagatcatgagtgtaaaaaactggaaatccaaagcctgaaatggctgccactcagaactgtttaagacattattgattccaagacagga
gaaacagcaagtaaacgatggaaaggtgccaaaaatagtgaaacagctgcaaaaggtgcatgtatgaaattaaagatgcatgctgatggcgataagtcattaccacagacagaaagaatttacttccaggttttcttaccataaaggagcaagagaagacaaaccaatgtcttttgcac
cgaatggagcattggaagggccatagacttgcgcttctagccaggcttaaaatgacaataacaaattacagctaagaaattaaggctgtgtcacattactcaggagaagccttaccctggatcatacttggaaacctggattgctaaggaggattgtcctttatataatgtgggaaat
aatcttggataatctaatgatgaagaacaattctgtaaaatgtgaatcttacttggaaatag