

SPC25-SPC25, NDC80 kinetochore complex component, homolog (*S. cerevisiae*) Gene

Reference: PTXBC022255

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SPC25, NDC80 kinetochore complex component; SPC25, NDC80 kinetochore complex component, homolog; kinetochore protein Spc25; AD024; SPBC25; hSpc25; 2600017H08Rik; spindle pole body component 25 homolog		

For research use only.

GENE DESCRIPTION	SPC25, NDC80 kinetochore complex component, homolog (<i>S. cerevisiae</i>)
GENE ID	57405
ORIGIN SPECIES	Human
ACCESSIONS	BC022255
VECTOR	pDONR223
NCBI SYMBOL	SPC25
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtagaggacgaactggcacttttcgataaaagcataaatgaattttggaataaattcaaaagtagcggacacctctgtcagatggcgggactaagagatactacaaggattccatcaaaagcatttgcagaaaagctgtctgtgaaattaaaggagaagaacgaatggttgagatgtttct
ggaatafcaaaatcagatcagcaggcaaaataagctcattcaagaaaaaaggataactgttaaaattgattgctgaagtaaaaggcaaaagcagggaattggaagtactgactgcaaatatccaggatcttaaggagaatattctaggaagaaggaaactatttctactctaataaagcg
aatgcagagaggttgaaaaggctgcagaaatctgcagactgtataaagatcgacttgactagaattcgaaaaatttatgggtgagaattgcagtttatttctactaatattgacccctaagaatcctgagagccatttatgttttcttactatcctcaatgaagcaagggaactatgaagtgcagat
agtccccctcatcttgaggccctagcagaattcaagagaatgaaggagaaccaacaattttcagctttcttgccaatgttcgaaagcttttactgccacggtttataattaa