

LGALS14-lectin, galactoside-binding, soluble, 14 Gene

Reference: PTXBC022257

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CLC2; PPL13; placental protein 13-like; Charcot-Leyden crystal protein 2; gal-14; lectin, galactoside-binding, soluble, 14; placental protein 13-like protein; galectin 14		

For research use only.

GENE DESCRIPTION	lectin, galactoside-binding, soluble, 14
GENE ID	56891
ORIGIN SPECIES	Human
ACCESSIONS	BC022257
VECTOR	pDONR223
NCBI SYMBOL	LGALS14
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcatcactaccgtaccatacacactgcctgtttccttgccctgttggttcgtgcgtgataatcacaggacacgcacatcctcactttgtcaaggacccacagctggaggtgaatttctacactgggatggatgagactcagatattgctttccaattccgactgcactttggcatcctgcaatc
atgaacagtcgtgtgttgcatatggagatagaggagaaatgctactatttacctttgaagatggcaaacatttgagctgtgcactatgtgcgtcacaaggaatacaagtaatgtaaatggccaacgcatttacaactttgccatcgattcccgccagcatctgtgaagatgctgcaa
gtcttgagagatattcctgaccagagtgcttatcagcgattga