

CD300C-CD300c molecule Gene

Reference: PTXBC022279

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD300c molecule; CD300c antigen; CLM-6; CMRF-35; CMRF-35A; CMRF35; CMRF35-A1; CMRF35A; CMRF35A1; IGSF16; LIR; CMRF35-like molecule 6; CD300 antigen-like family member C; CMRF35 antigen; CMRF35 leukocyte immunoglobulin-like receptor; CMRF35A leukocyte immu		

For research use only.

GENE DESCRIPTION	CD300c molecule
GENE ID	10871
ORIGIN SPECIES	Human
ACCESSIONS	BC022279
VECTOR	pDONR223
NCBI SYMBOL	CD300C
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgactgccaggccctggcctcgtggcgtcttcagctctgctcctctgtgtccaggtattttctctgagccacccatgacctggcgggccccgtgggggatccctgagtgctgctgctatgagaaggacacaggacctcaacaattctggtcagaccaccac
agattctccgatgtgacaagattgtgagaccaaaagggtcagcaggggaaagggaatggccgagtgctccatcaggacagctctgcaaacctcagcttcacagtgacctggagaatctcacaggaggacgcaggcacctactgtgtgtgggtggataccctggtccgagacttt
catgatcccatgtcaggttgaggtgtccgtgttccggccgggacgacacacgctccagccccagagctccatgggcacctcaggtctccacgaagctgcccgtgcacacctggccagcgtgacgagaaggacagccccgaaccagccacacctggtccctgttca
gcaatgtccgtctctgctcctgtgagctgcccctgctcctgagcatgtgggtgccgtcctctgggtgaacagacctcagagaagcttagaagcaggcagaattggccaagggtgagaaccagtag