

CFHR2-complement factor H-related 2 Gene

Reference: PTXBC022283

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CFHL2; FHR2; HFL3; complement factor H-related protein 2; DDESK59; FHR-2; H factor (complement)-like 3; factor H-related gene 2; h factor-like 3; h factor-like protein 2; complement factor H related 2		

For research use only.

GENE DESCRIPTION	complement factor H-related 2
GENE ID	3080
ORIGIN SPECIES	Human
ACCESSIONS	BC022283
VECTOR	pDONR223
NCBI SYMBOL	CFHR2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgctcctggtcagtgtaattctaattcagcgatattctgtgtgggggagaagcaatgttctgtgattttccaaaaataaaccatggaattctatatgatgaagaaaaataagccattttcccaagttctacaggggaagttttctattactcctgtgaatataatttgtgtccttcaaaat
cccttttggactcgataacgtgcgcagaagaaggatggtcaccaacaccaaaggtgtcagactgtgtttctttcttttgggaaaaatggtcattctgaattctcaggacaaacacatctggaaggtgatactgtacaattatttgaacacaggatagacttcaaaacaatgagaacaacat
ttcatgtgtagaacgggctgtgtccactctcccaatgcaggtccactatttctgcagaaaaatgtgggccccctccactattgacaatggagacattacttctcgtgtgtcagtatatgtccaggttcatcagttgagtagcaggtgccagaactgtatcaactgagggtaacaacaaa
taacatgtagaacggacaatggtcagaaccaccaaattgcttagatcctgtgtatatacacaagaaattatgaaaaataatacataaaataaagtggacaaccaacaaaagctttatcaagaacaggtgacatagttgaattgtttgtaaatctggatatcatcaacaaaatctcattca
tttcgagcaatgtgcagaatgggaaactggtatattccagttgtgaagaaaatag