



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

RAMP3-receptor (G protein-coupled) activity modifying protein 3 Gene

Reference: PTXBC022304

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

receptor activity-modifying protein 3; CRLR activity-modifying protein 3; calcitonin receptor-like receptor activity modifying protein 3; receptor (G protein-coupled) activity modifying protein 3; receptor (calcitonin) activity modifying protein 3; recept

For research use only.

GENE DESCRIPTION	receptor (G protein-coupled) activity modifying protein 3
GENE ID	10268
ORIGIN SPECIES	Human
ACCESSIONS	BC022304
VECTOR	pDONR223
NCBI SYMBOL	RAMP3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagactggagcgctgcggcgcccgcaactctccgtgtgctgctgctctgggggtgccagagcaggcggtcgcaacgagacagcagctgttggagagggctgccccgtgtggggaaggcttgcagacatgatggcgcaaggctgacgtctggaagtgtgcaactctgcgaggtctatcgtgtgactatgagagtttcaacgcacccagagatgagggcgaactgtgtggcgtctactgtccctcccaacccccctggccagggcttcatcacccgtcatccacagcgagttgtcttccaaactgcacgttgagacagggctcacttggagacccccagacagaggttctcat