

SERPINC1-serpin peptidase inhibitor, clade C (antithrombin), member 1 Gene

Reference: PTXBC022309

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AT3; AT3D; ATIII; THPH7; antithrombin-III; serine (or cysteine) proteinase inhibitor, clade C (antithrombin), member 1; serpin peptidase inhibitor clade C member 1; serpin peptidase inhibitor, clade C (antithrombin), member 1; serpin family C member 1		

For research use only.

GENE DESCRIPTION	serpin peptidase inhibitor, clade C (antithrombin), member 1
GENE ID	462
ORIGIN SPECIES	Human
ACCESSIONS	BC022309
VECTOR	pDONR223
NCBI SYMBOL	SERPINC1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtattccaatgtgatagggaactgtacactctggaaaaaggaaggttatctttgtccttgctgctcattggtcttgggactgctgacgtgacgggagccctgtggacatctgcacagccaagccgcccgaattcccatgaatccatgtgcattaccgtccccggagaagaaggaagc
aactgaggatgagggtcagaaacagaagatcccgaggccaccaaccggcgtgtctgggaactgtccaaggccaattcccgctttgtaccatttctatcagcactggcagatccaagaatgacaatgataacatttctgtcaccctgagtatctccacggctttgtctatgaccaag
ctgggtgctgtaatgacacctccagcaactgatggaggttaattgaacaccatactgagaaaaacatctgacagatccacttctcttgcacaaactgaactgccgactctatgaaaagccaacaaatctccaagttagtagatagccaatcgcttttggagacaaatcccttacctt
caatgacctctatgtctcagatgcatccataaggcatttctgaggttaaatgaagaaggcagtgaaagcagctgcaagtaccgctgtgtgattgctgcccgttcgctaaacccaacaggggtgacttcaaggccaacaggcccttctgtgtttataagagaagttcctctgaacactattatct
tcatgggcagagtagccaacctgtgttaagtaa