

C16orf63-chromosome 16 open reading frame 63 Gene

Reference: PTXBC022321

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	lisH domain-containing protein C16orf63; C16orf63; FOR20; PHSECRG2; lisH domain-containing protein FOPNL; FOP-related protein of 20 kDa; pluripotent embryonic stem cell-related protein; FGFR1OP N-terminal like		

For research use only.

GENE DESCRIPTION	chromosome 16 open reading frame 63
GENE ID	123811
ORIGIN SPECIES	Human
ACCESSIONS	BC022321
VECTOR	pDONR223
NCBI SYMBOL	C16orf63
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgactgtggcagagttgaaggctgttttaaggacaccttggaaggggggtattagggcatttaaaagcaaggatccgagctgaagtttcaatgccctagatgacgtgaaccccgaccatcattgtctcatgaaaccttctaataatgaattaattcgagagtatttagaa
ttcaacaaatataagatagcagcatctgtctcatagcagaatctgtcaacctgtattccgttggacagacagtttctcatcgaactaaatgcatttgaagaatcaaaggataataacacctctttatattgggatttttagccatttcttgcgtggaaactaaggatggcatccagaatgcatt
tctgaaaggccttcacttcagccttcagaccaagcttggcagacaacctagtagaagaagccaatggatgaccacctaagaaggagggaacagaaaagtactaacaattgaagatcttcatgtttctcaggcagtcacagataa