

PHOSPHO2-phosphatase, orphan 2 Gene

Reference: PTXBC022324

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	pyridoxal phosphate phosphatase PHOSPHO2; phosphatase, orphan 2		
For research use only.			

GENE DESCRIPTION	phosphatase, orphan 2
GENE ID	493911
ORIGIN SPECIES	Human
ACCESSIONS	BC022324
VECTOR	pDONR223
NCBI SYMBOL	PHOSPHO2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaaatttgcagttttgactttgacaatacatagatgacaatagtgacacttggattgtacaatgtgtcccaacaaaagcttctattgaactacgtgattcttatcgaaaaggatttggacagaatttatggcagagctttaaagtattgggagataagggtgaagagaacatga
aatgaaaagagcagtgacatcattgcttctactccagggatgggtgaactctcaactttataagaagaataaggataaattgactgcattattttcagatcaaatcggcttcatagattgggtttagaagctgccagtttcatgacatattgataaagtgtttacaatccagcagctttt
aatagcaatgggtcatctcactgttgaaaattatcatactcattcttgcataatagatgcccaagaatcttgcacaaaaggtagttttgatagaattttagataaacagttacaacagggagtgattatacacaattgtttatatgtgtgagtggtgaaatgatgtctgtccagtcaccttttaagaat
gatgatgttgcacgggaaaggatatacttacagaaaactcttccagaatgtctcaaatcttgagcctatggaattattctgttagtttggctcaggtgttgatataatttctcattacaatttctaataaaggattaa