

PHC3-polyhomeotic homolog 3 (Drosophila) Gene

Reference: PTXBC022325

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	EDR3; HPH3; polyhomeotic-like protein 3; early development regulator 3; early development regulatory protein 3; homolog of polyhomeotic 3; polyhomeotic like 3; polyhomeotic homolog 3		

For research use only.

GENE DESCRIPTION	polyhomeotic homolog 3 (Drosophila)
GENE ID	80012
ORIGIN SPECIES	Human
ACCESSIONS	BC022325
VECTOR	pDONR223
NCBI SYMBOL	PHC3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggaagcgggaatttaaggaccatagtagcagctatggatactgaaccaaacccgggaacatcttctgtgtcaacaacaaccagcagtagcaccaccaccaccatcaccacttctctctcgaatgcagcagccacagatctctgtctacagtgggtcagaccgacatgctgtacaggt
aattcaacaggcattgcatcgccccccagctcagctgtcagtagcttcagcaaatgtatgcagcccaacaacagcacttgatgctgcatctcagcagcagcatttaagcagctccagcttcagagccttgcgtgttcaggcaagttgtccagtggaagaccatctacat
ctccacaggaagtgtagcacagcagtagtcaagtagtcccaacgtctgtaagttccctaaatttttttctcctggatttaaaatttttaa