

PLAC1-placenta-specific 1 Gene

Reference: PTXBC022335

| | | | |
|----------------|--|---------------------|----------------|
| Product type: | DNA | Storage Condition: | -20°C |
| Size/Quantity: | 2 ug | Delivery Condition: | +4°C, Blue Ice |
| Synonyms | CT92; OOSP2L; placenta-specific protein 1; cancer/testis antigen 92; placenta specific 1 | | |

For research use only.

| | |
|-------------------------|---|
| GENE DESCRIPTION | placenta-specific 1 |
| GENE ID | 10761 |
| ORIGIN SPECIES | Human |
| ACCESSIONS | BC022335 |
| VECTOR | pDONR223 |
| NCBI SYMBOL | PLAC1 |
| LENGTH | ProteoGenix bp |
| SEQUENCE PRIMERS | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

ORF (Only for reference, product provided without STOP codon)

atgaaagttttaagttcataggactgatcctcctacctctgcgtttcagccggttcaggacaaagccaatgactgtgctgtcctcatagactggttcacagtgccaccctcatgctaaacaacgatgtgtgtacactttcatgaactacactgggctgggtgccccca
aaccatgtcagccacacgcctaccagttcacctaccggtgtactgaatgtggcatcaggccaaagctgtctctcaggacatggtttactacagcactgagatacactacttctaaggcagccatctaaagttgtgatcccagtgcatgtgctgcccccaaaagtccccatggctcac
caagccctgctccatgagagtagccagcaagagcagggccacagcccagaagatgagaaatgctcagaggtgtcagctgtcacagtcagtcacaaaggcccaactgcgattgtccacctgtgtctcagtggaagaagagcataccaggtccctgtcaccaagcagggctcagg
aggctcaacctctgcagccatcactttctgatattctgaggattgtctctcacacagatgatgattgggtccatgta