

C18orf32-chromosome 18 open reading frame 32 Gene

Reference: PTXBC022357

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UPF0729 protein C18orf32; chromosome 18 open reading frame 32		
For research use only.			

GENE DESCRIPTION	chromosome 18 open reading frame 32
GENE ID	497661
ORIGIN SPECIES	Human
ACCESSIONS	BC022357
VECTOR	pDONR223
NCBI SYMBOL	C18orf32
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtgtgcattcctgtatcgtcattccagttctgctctggatctacaaaaattcctggagccatatatataccctctggtttccccccttgtagtcgtatgacctagaaagcaatacaagaatccaatgatacaacaaggcaagtaactttaagggtgcagacatgaatggattacca
acaaaaggaccaacagaaatctgtgataaaaagaagactaa