



Reference: PTXBC022367

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CRISTIN1; PWTSR; THSD2; R-spondin-3; R-spondin 3 homolog; protein with TSP type-1 repeat; roof plate-specific spondin-3; thrombospondin type-1 domain-containing protein 2; thrombospondin, type I, domain containing 2; R-spondin 3		

For research use only.

GENE DESCRIPTION	R-spondin 3 homolog (<i>Xenopus laevis</i>)
GENE ID	84870
ORIGIN SPECIES	Human
ACCESSIONS	BC022367
VECTOR	pDONR223
NCBI SYMBOL	RSPO3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]