

BTF3L4-basic transcription factor 3-like 4 Gene

Reference: PTXBC022371

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	transcription factor BTF3 homolog 4; basic transcription factor 3 like 4		

For research use only.

GENE DESCRIPTION	basic transcription factor 3-like 4
GENE ID	91408
ORIGIN SPECIES	Human
ACCESSIONS	BC022371
VECTOR	pDONR223
NCBI SYMBOL	BTF3L4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaatcaagaaaagttagcacaacttcaggctcaggtccgataggggcaagggtacagctcgcagaaaagaaggtggtacatagaacacccacagctgatgacaaaaagcttcagagttctctaaaaaaactggctgtgaataatagctgtattgaagagtgaaatgatta
aagatgatgggacagttattcattcaacaatcccaagtccaagctccctttctgtaatacctttgcaattactggtcatgcagaagccaacaacacagaaatgcttctggaatattaagtcagcttggctgacagttacaagccttaggaagtgctgaaacagttcccacgcaa
gtcttgacagtaaaacacaaaaccagaagacattgatgaggaagatgatgatgtccagatctgtagaaaatttgatgagcatcaagaatgaagctaactaa