

PHF19-PHD finger protein 19 Gene

Reference: PTXBC022374

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MTF2L1; PCL3; TDRD19B; PHD finger protein 19; polycomb like 3; polycomb-like protein 3; tudor domain containing 19B		

For research use only.

GENE DESCRIPTION	PHD finger protein 19
GENE ID	26147
ORIGIN SPECIES	Human
ACCESSIONS	BC022374
VECTOR	pDONR223
NCBI SYMBOL	PHF19
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagaatcgagctctggatccaggactcgggactcctatgggtgccaccagccactccccacaagggggccctggcggaaggtcaagaacaactcaaaagacttgatgtccaaactgacggaggccagtagtgcgtgccggtggacagatggcctgtactacctgggaaga
tcaagagggtcagcagctctaaagcaagctgcctcgtgacttgcgaagataatccaaatactgggtcctatggaaggacatacagcatgccgggtgttccaggagagagcccaagtgaacatctgcctagggaagacatcaggccgctgaatgagatcctcatctcgggaagtgtg
gcctgggttaccaccagcagtgccacatccccatagcgggcagtgctgaccagccctgtcacaccttgggtctgccgacgtgcattctgcactgctgtgcgggtgagccttccatctccccagtcctgcctcctcctcagtggggcagaccagagactccatcacag
agtctgagctccaagcagaaggccacactgggttggagacagatagcgcctctgcactgtccttggccaggattgttag