

C4orf39-chromosome 4 open reading frame 39 Gene

Reference: PTXBC022381

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C4orf39; C4orf39; protein FAM218A; family with sequence similarity 218 member A		

For research use only.

GENE DESCRIPTION	chromosome 4 open reading frame 39
GENE ID	152756
ORIGIN SPECIES	Human
ACCESSIONS	BC022381
VECTOR	pDONR223
NCBI SYMBOL	C4orf39
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaaggatgcggtgctggcgccggcagctgtctcttccccggaccagcgccctggagagccagccctgcagggtgggtggcgagccaaactcgttctgggtgcagggtcctgggtctcctaacagacctatacgtgacggcgccgccatggcagtgctctttgtctcagacatccaggagaccacattcgtccaacagcgggtgctccaccaatcctgggagagcgaatcgttttctcgcgtgccctgtcagccgtcatgtgtcccagagaggaaatttagtgagcattccgggtgtcagccaccgaaatfgccaggccactccaagtcagaaggaccaccaggaaaagtcaggaagagaaccacatcagggtccagcctcttttggacaaggactagagggtttgggtctgcagtggtggtcctcccgctgggtcacctgtcctatag