

MRPL30-mitochondrial ribosomal protein L30 Gene

Reference: PTXBC022391

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	L28MT; L30MT; MRP-L28; MRP-L30; MRPL28; MRPL28M; RPML28; 39S ribosomal protein L30, mitochondrial; 39S ribosomal protein L28, mitochondrial; mitochondrial ribosomal protein L30		

For research use only.

GENE DESCRIPTION	mitochondrial ribosomal protein L30
GENE ID	51263
ORIGIN SPECIES	Human
ACCESSIONS	BC022391
VECTOR	pDONR223
NCBI SYMBOL	MRPL30
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctgggatttgccttagttagtcaatggccccagcagactacagaccgtgacaaaagggtggaggtctcttattgtacagattggattcgtcacaattcaccagatcaagaattccagaaaaagtggttcaggcctcacctgaagatcatgaaaaatacgggtgggatccacaga
accctcataaactgcattgttaccagaataaaaagtacaagaagacgtccatattgggaaaaagataataaagatgcttgattagaaaaagcacatcccctcaagttcacaagaatcccttcagtgaaatgcaaaattgaaagtagttaagcatttgataagaatcaagcccttgaagtt
gccacaaggacttccagcagaggagaacatgtctaacacgtgcctcaaaagcactggggagtttagttagtcagtggtgcatctgaaacctgtggagcagaaagcacatgagtcctaa