

RICH2-Rho-type GTPase-activating protein RICH2 Gene

Reference: PTXBC022452

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	rho GTPase-activating protein RICH2; Rho-type GTPase-activating protein RICH2; RICH2; NPC-A-10; rho GTPase-activating protein 44; RhoGAP interacting with CIP4 homologs protein 2; Rho GTPase activating protein 44		

For research use only.

GENE DESCRIPTION	Rho-type GTPase-activating protein RICH2
GENE ID	9912
ORIGIN SPECIES	Human
ACCESSIONS	BC022452
VECTOR	pDONR223
NCBI SYMBOL	RICH2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcagaccagtcgctggccagccgtcccccagtcagcctgtccccacccccccagcagccccgtaccctatggactgagctaccctcaggggtactccttggcctcgggccagctctcccccagctgcagctcctccctggcctcctctgtctttacaagcactttgagcaaatc
gcggcccaactcctaagccgcgacagagacctactctgccgctcctcagcctccacagtaaacctctcggcctctagtcacagtcacaggagggcccatgctagatggcatgtccctggggaaagcatgtctacagctgtgtga