

COQ10B-coenzyme Q10 homolog B (*S. cerevisiae*) Gene

Reference: PTXBC022453

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	coenzyme Q-binding protein COQ10 homolog B, mitochondrial; coenzyme Q10 homolog B; coenzyme Q10B		

For research use only.

GENE DESCRIPTION	coenzyme Q10 homolog B (<i>S. cerevisiae</i>)
GENE ID	80219
ORIGIN SPECIES	Human
ACCESSIONS	BC022453
VECTOR	pDONR223
NCBI SYMBOL	COQ10B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcagctcggactggatcacggccttgagaaggtagtctcgggatccgtccgaagtcggcagacagcggccggcgaggcggccgtgcggaatggcagatatttagcttcctgtgtgtatctgatgagcagaactcttcactacatacctcaatttgcctaaggagatatg
cacgaactttctcaaaatcactgcaccattaataacaaaaggaaagaatttcagagagaagaatttttagatattcaatgcaggaaatgtatgatgtatcgggagtgaggattacaagcattttgtccttggtgcaaaaaatcagatgttatataaagagatcggatattgtaaaaca
agattagaaatggattccacctgtgttgagcgatatacatcagtagtaacctgtgtgaaacctcaatttagtaaaaggcatctgtactgatgggagactttcaatcattggagactattggcgttttagcccggtcttcctggctaccgaagaactgtacctggattttcaattctttgaattt
cgatcacttctacattccagctgccacactctttttgatgaagttgtgaagcagatgtagctgcctttgaaagaagagcatgtaagctgtatgtccagaacaataatacctcggagttaatgctcatgaagtcacatcacataa