

## FEZF2-FEZ family zinc finger 2 Gene

Reference: PTXBC022464

|                |                                                                                                                                                                                                                         |                     |                |
|----------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                                                     | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                                                    | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | FEZ; FEZL; FKSG36; TOF; ZFP312; ZNF312; fez family zinc finger protein 2; forebrain embryonic zinc finger-like protein 2; testis tissue sperm-binding protein Li 80P; zinc finger protein 312; FEZ family zinc finger 2 |                     |                |

**For research use only.**

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | FEZ family zinc finger 2                                                                         |
| <b>GENE ID</b>          | 55079                                                                                            |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC022464                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | FEZF2                                                                                            |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon)

atggcaagctcggttccctggagaccatgggtgccccggcctgcccgcgcgcggagcgtgcggcgccacttcaagacactggcctttccatcgagcgcacatggccaagacgtcgaggccccgtgcgcctttgagccccggcctggagcgtagaggcgagcgagcca  
gggcaagaaactgctcaacctctgctcgccgtgcccgtgatgacccctccagccctaggctacgaggtgcccgtcaagacactgctcagttactcggagctctggaagcagcctccggcgggcgggcgggcgaggcgggcggtggcgggcgggcgggggg  
ggccccagtggtgcggcgccagcggcgtgtgtgcaaaaccaactgtggcgtgtgtgtgcaaggccgagctggcctggcgccgtccgcgtgcccggcgagggtcatcaagccgaggtcatcaaccagcgtgtgggctgacggccagcggcgtcgtctactactcaactacgtg  
actcgaccgcgtaccgccgtgtgagctcctcagcggccacactctcccgctgtgcccctcctcaatgcgcagggcccccgcccgccctggtcgtcaccaccaagctctttctgtggagaacgccaaagctggccggcctggtcgggacaagttccccaccggcctccatccccataagg  
agcgttgcggcgccgtggagcaggtactgaaggaactggccctgactgcccagcgcggagcgtcaaggccacagcaagctccaggaggtccgcagatggcagcccaaaacttccctgcgaggtgtgtggcaaggtgttaacgtcactataatctacccgc  
cacatgccgtccacaccggagcagaccgttctgtgcaagctgtgcgcaaggcttccgagcgccagcagcgtctgcagccacaaaattatccacaccaggaaaagccacataaatgcaaccagtcggcgaagcgttcaaccgcagctccacgtcaacacgcataatccga  
tccacgcgggtacaaagcccttcgtcgtgcaattttgcggcaagggtttccacaaaaagggaactacaagaaccacaagctgacccacagcggcgagaagcagtagacaatgtaccatctgcaacaaggccttccaccaggtctacaacctaaccttcataatgcacaccacaacgaca  
agaagcctttcacgtgcgccacttgcggcaagggtttgcagaaacttgactaaagaacatgtgcgcaaaactccacgacagcgtggccctgctgccccctccgcaaggacctgactaggacagtgacagctga