

## INSIG2-insulin induced gene 2 Gene

Reference: PTXBC022475

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	INSIG2 membrane protein; INSIG-2; insulin-induced gene 2 protein; insulin induced protein 2; insulin induced gene 2		

**For research use only.**

---

<b>GENE DESCRIPTION</b>	insulin induced gene 2
<b>GENE ID</b>	51141
<b>ORIGIN SPECIES</b>	Human
<b>ACCESSIONS</b>	BC022475
<b>VECTOR</b>	pDONR223
<b>NCBI SYMBOL</b>	INSIG2
<b>LENGTH</b>	ProteoGenix bp
<b>SEQUENCE PRIMERS</b>	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

### ORF (Only for reference, product provided without STOP codon)

atggcagaaggagagacagagtcacctggcccaaaaagtgtgcccataatfcatctgactagccagagtgtgaactgatgattcgaggagtgtgtattttattggagtattcttcattagtgttaattfactcagattcagagaatgtgacgctttccacctgatgtgattg  
caagcatctttctctgcatgggtgacccccatgctgtggcagggctcagctgtgattgggtattatacccctgcattgacagacatclaggagaaccacataaattaaaagagagtgtccagtgtaatgcgggtgtgtagcagctttgtgtataaatcatgccagtctaaagtgatt  
tcgataacaacatacagttgtctcacactggctgcaclatccattggactgtgtggactttgalagatctagaagtgtttggcctggagtaggaattgccctttggcaactgtgtcactcaactgctagtataatagggtttaccataatcatctccagattcctctatgtctcttg  
gttaccatgtatattttgtggaggcataacaatgggaacattgtcgacaactggcaatgtacgaatgtaagtattcgcagaaaaatctcatcaggaatga