

FABP6-fatty acid binding protein 6, ileal Gene

Reference: PTXBC022489

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	I-15P; I-BABP; I-BALB; I-BAP; ILBP; ILBP3; ILLBP; gastrotropin; fatty acid binding protein 6, ileal; ileal bile acid binding protein; ileal lipid-binding protein; ileal lipid-binding protein; intestinal 15 kDa protein; fatty acid binding protein 6		

For research use only.

GENE DESCRIPTION	fatty acid binding protein 6, ileal
GENE ID	2172
ORIGIN SPECIES	Human
ACCESSIONS	BC022489
VECTOR	pDONR223
NCBI SYMBOL	FABP6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcttcaccggcaagttcgagatggagagtgagaagaattatgatgagttcatgaagctccttgggatctccagcagatgaatcgaaggcccaactcaagatcgtcacggaggtgcagcaggatgggcaggactcacttggtccagcactactacggggccacacatga
ccaacaagttcactgttggcaaggaaagcaacatacagacaatggggggcaagacgttcaaggccactgtgcagatggagggcgggaagctgtgtgtgaatttcccaactatcaccagacctcagagatcgtgggtgacaagctgtgtggaggtctccacatcggagcgtgaccta
tgagcgtgagcaagagactggcctaa