

CD200-CD200 molecule Gene

Reference: PTXBC022522

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD200 molecule; CD200 antigen; MOX1; MOX2; MRC; OX-2; OX-2 membrane glycoprotein; antigen identified by monoclonal antibody MRC OX-2		

For research use only.

GENE DESCRIPTION	CD200 molecule
GENE ID	4345
ORIGIN SPECIES	Human
ACCESSIONS	BC022522
VECTOR	pDONR223
NCBI SYMBOL	CD200
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagaggtggtgatcaggatgccctctgtcatctgtctacacagcctgtttgggtcatggcagcagtggtgtgtgtcacagcacaagtgcagtggtagccaggatgaaagagagcagctgtacacacctgtctcttaaaatgtctctgcataatgccaggagccctcatt
gtgacatggcagaaaaagaaagctgttaagcccagaaaacatgtgtcaccctcagcgagaacatgggtgtgtatccagcctgacctataaggacaagataaacattaccagctgggactccaaaactcaaccatcaccctctggaatatcaccctggagatgaagggtgttacatgtgtc
tcttcaatacctttgttttgggaagatctcaggaacggcctgctcaccgtctatgtacagccatagatcccttactacaaaatctctgaagaccacctaataatcactgtctgtccactgccgcccagccccatgggtcttctggaaggtccctcgggtcaggattgaaaatagtagt
gactgtgtcaccacaaatgggaccacgtgtgtaccagatcctccatcaaaagaccctaagaatcagtggtgggaaggaggtgatctgccaggtgtgtcacctgggactgtgaccgacttaagcaaacctgcaacaaggctattgttttcagttccgtattgctaagcattgttccc
tggttaattctctgtcctaattctcaattctgtactgtgaaacgtcaccggaatcagagaccgagagccctaa