

ZCRB1-zinc finger CCHC-type and RNA binding motif 1 Gene

Reference: PTXBC022543

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MADP-1; MADP1; RBM36; SNRNP31; ZCCHC19; zinc finger CCHC-type and RNA-binding motif-containing protein 1; U11/U12 small nuclear ribonucleoprotein 31 kDa protein; U11/U12 snRNP 31 kDa protein; U11/U12 snRNP 31K; U11/U12-31K; zinc finger CCHC-type and RNA b		

For research use only.

GENE DESCRIPTION	zinc finger CCHC-type and RNA binding motif 1
GENE ID	85437
ORIGIN SPECIES	Human
ACCESSIONS	BC022543
VECTOR	pDONR223
NCBI SYMBOL	ZCRB1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

```
atgagtgtggattggctccaagtaagacacagtgatgtalccaactgacctttccctgacaacaatgactgtaccggatatttccaagatggcaagttgtaagggtaccatcatgaaagataaagataaccaggagagtaaggggtgcattatattttggataaagactctg
cacaaaactgaccaggccaataaacaacaacagttatftgtagagtataaaagcaagcattgctattgacaatggaagagcagctgagttcatccgaagcgaaactacttgataaatctaagtgttatgaatgtgggaaagtgacactaaagttatgcctgcccgaataatgctc
ggagaactgtagcctcaaaagaagaagaaaaaaagaaaaaaagaaagctcctgaaccagaagaagaattgaggaagtagaagaagtgaaatgaaagggaggatcctgctctgacagcctcagtcaggccatagcattccagcaagccaaatgaaagaacaaaaaa
atggaaccagttcaggatccctcaacatcagatgattcaagacgcccaaggataaagaaaacacatatttcagtgatgaggaagaactagtgattaa
```