

ABHD13-abhydrolase domain containing 13 Gene

Reference: PTXBC022566

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein ABHD13; BEM46L1; C13orf6; bA153I24.2; abhydrolase domain-containing protein 13; alpha/beta hydrolase domain-containing protein 13; abhydrolase domain containing 13		

For research use only.

GENE DESCRIPTION	abhydrolase domain containing 13
GENE ID	84945
ORIGIN SPECIES	Human
ACCESSIONS	BC022566
VECTOR	pDONR223
NCBI SYMBOL	ABHD13
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttggttaacctcaagtttaaccttttctggttgattatcgaggataggaaaaagtgaaggagaagcaagtgaagaaggactctacttagattctgaagctgtgttagactacgtgatgactagacctgaccttgataaaacaaaaattttcttttggccgttccttggtggagcagtggtct
attcatttggcttctgaaaaatcacataggatttcagccattatggtggagaacacatttttaagcatalccacataggccagcactttatttattcttccgatgcgttaccttctttatggtgtacaaaaataaattttgtctacagaaaaatctctcagtgtagaatgccttcaacttttcatctctgg
actctcagatcaatttaattccaccagtaattgatgaacaactttatgaactctcccatctcggactaagagattagccattttccagatgggactcacaatgacacatggcagtgccaaggctatttcaactgcacttgaacagttcatcaagaagtcgtaagagccattctctgaagaaatg
gcaaaaactcatctaatgttaacaattatataa