

C5orf40-chromosome 5 open reading frame 40 Gene

Reference: PTXBC022570

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	fibronectin type-III domain-containing protein C5orf40; C5orf40; fibronectin type III domain-containing protein 9; fibronectin type III domain containing 9		

For research use only.

GENE DESCRIPTION	chromosome 5 open reading frame 40
GENE ID	408263
ORIGIN SPECIES	Human
ACCESSIONS	BC022570
VECTOR	pDONR223
NCBI SYMBOL	C5orf40
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaacatcagagtgagggaacattcttatacaggagccatcctctgtcgtcctcggagccctgcctggaggaactattaccatattatgtacaggccaactggaacagcatcttctctggtatcttctgctacagctccacaacgaggagaaggtgcctcgaacgatcagctcgtgg
tgcctggaacatcttgcccttccactctactctctgtgcatcagctgtaagaaggctgcttcccttacaggcactactgcacatgtccacacccctggataagagtcctgctgctggaagctccctgtagacccccagatctcccttgggtgctgatggccattctgctgctgctt
cacagccgtcttgccctcattctgctccagttctgtgtgctgctgctgcatgagccgcgatggcttctacaggcctggccatggaggagccaatgggtggtgagatggccagaggagggcccgatcttggctcagaggagggaagacctgcaggggtccccctggtggaatgc
cacgcaagaactccagagatggagctgaactggatccggaagccaaccaggatgccccgatgcgggtgccttacagaggggggtggtgaccaccgcctatactgcctcattgtgggaatga