



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

CGB-chorionic gonadotropin, beta polypeptide Gene

Reference: PTXBC022796

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CGB; CGB5; CGB7; CGB8; hCGB; choriogonadotropin subunit beta 3; CG-beta; choriogonadotropin subunit beta; chorionic gonadotrophin chain beta; chorionic gonadotropin beta 3 subunit; chorionic gonadotropin beta chain; chorionic gonadotropin beta subunit; ch		

For research use only.

GENE DESCRIPTION	chorionic gonadotropin, beta polypeptide
GENE ID	1082
ORIGIN SPECIES	Human
ACCESSIONS	BC022796
VECTOR	pDONR223
NCBI SYMBOL	CGB
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr