



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

CCDC28B-coiled-coil domain containing 28B Gene

Reference: PTXBC022848

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	coiled-coil domain-containing protein 28B; coiled-coil domain containing 28B		
For research use only.			

GENE DESCRIPTION	coiled-coil domain containing 28B
GENE ID	79140
ORIGIN SPECIES	Human
ACCESSIONS	BC022848
VECTOR	pDONR223
NCBI SYMBOL	CCDC28B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]