

CYP4A11-cytochrome P450, family 4, subfamily A, polypeptide 11 Gene

Reference: PTXBC022851

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CP4Y; CYP4A2; CYP4AII; cytochrome P450 4A11; 20-HETE synthase; 20-hydroxyeicosatetraenoic acid synthase; CYP1VA11; P450HL-omega; alkane-1 monooxygenase; cytochrome P-450HK-omega; cytochrome P450, family 4, subfamily A, polypeptide 11; cytochrome P450, sub		

For research use only.

GENE DESCRIPTION	cytochrome P450, family 4, subfamily A, polypeptide 11
GENE ID	1579
ORIGIN SPECIES	Human
ACCESSIONS	BC022851
VECTOR	pDONR223
NCBI SYMBOL	CYP4A11
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagtgtctctgtctgagccccagcagactcctgggtgatgtctctggaatctccaagcgccctcctgtcattctgtctctgtctgatcaaggcagttcagctctacctgcacaggcagtggtctctcaagccctccagcagttcccggtgccctccctccactggctctcgggca
catccaggagctccaacaggaccaggagctacaacggattcagaaatgggtggagacattcccaagtgctgtctcattggctatggggaggcgaagtctgtctcagctctatgacctgactatagagggtgattctggggagatcagacccgaaatcccatgggtctacagattcct
ggctccatggattgggtacgggtgtctctgtgaatgggcagacatggtccagcatcgacggatgtgaccccgagccttcactatgacatcctgaagccctatgtggggctcatggcagactctgtacgagtgtgctggacaaatgggaagagctcctggccaggattcccctctgga
ggtctttcagcagctctcctgatgacctggacaccatcatgaagtgtgccttcaggcattggcagagctcagcactccctcaccttcctga