

MS4A6A-membrane-spanning 4-domains, subfamily A, member 6A Gene

Reference: PTXBC022854

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MS4A6A-polymorph; 4SPAN3; 4SPAN3.2; CD20L3; CDA01; MS4A6; MST090; MSTP090; membrane-spanning 4-domains subfamily A member 6A; CD20 antigen-like 3; CD20-like precursor; HAIRB-iso; four-span transmembrane protein 3; four-span transmembrane protein 3.1; four-		

For research use only.

GENE DESCRIPTION	membrane-spanning 4-domains, subfamily A, member 6A
GENE ID	64231
ORIGIN SPECIES	Human
ACCESSIONS	BC022854
VECTOR	pDONR223
NCBI SYMBOL	MS4A6A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacatcacacacgtgtcccaatgagaccatcatagtgctccatcaaatgtcatcaacttctcccaagcagagaaacccgaaccaccaaccaggggcaggatagcctgaagaacatctacacgcagaaatcaaagtattgggactatccagatctgtgtgcatgatggtattgagc
ttggggatcattttggcatctgcttcctctccaaattttacccaagtgtactctacactgttgaactctgtctacccattcatagaccctttttttatcatctctgctctctatcaatgccacagagaaaagggtgaccaagctttgggcatagcagcctggttggagcattctgagtctct
gtctgccctgtgtgggtttcattatctgtctgtcaaacaggccacctaaatctgctcactgcagtgagttggacaaaaataataccaacaagaagtattgttcttactttatcatgattcatttataccacggactgtatagccaaagccagctgtggctggatccctctctgatgc
tgatttgactctgctggaattctgctagctgtctcactgtgtgtctgcgttggaacaggcttactctgacttccctggagtgacttttctgctcagttacattgtaattctgcatgtcctcaaaatgactcatgactgtggatagaagaactattgactcttaa