

XIRP2-xin actin-binding repeat containing 2 Gene

Reference: PTXBC022888

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CMYA3; xin actin-binding repeat-containing protein 2; beta-xin; cardiomyopathy associated 3; cardiomyopathy-associated protein 3; myomaxin; xeplin; xin repeat protein 2; xin actin binding repeat containing 2		

For research use only.

GENE DESCRIPTION	xin actin-binding repeat containing 2
GENE ID	129446
ORIGIN SPECIES	Human
ACCESSIONS	BC022888
VECTOR	pDONR223
NCBI SYMBOL	XIRP2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggatgcatttgagagtgcaaatgttgagtcgaagatgaaaacctcttcacatagctcagaagctggcaaatctgctgtgactcaagcatgccccaccaacctatgaggatgtcattgctggacatatttagatctctgattcacctaaagaagtaagaaaaatttcaaaagacgtg
gcaagagagtggaagagtttttaagcctgggatatgcaaccgcagatgcttctgcaactgagatgagaaccacctccaagaggaaatctgcatttataagtgaagctgctgctcgaagacaaggaaatatgtatactttgcaaaagacagtttatccaatggagtgccctagtgccagacaa
gcagaattttcataa