

MS4A8B-membrane-spanning 4-domains, subfamily A, member 8B Gene

Reference: PTXBC022895

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MS4A8B; 4SPAN4; CD20L5; MS4A4; membrane-spanning 4-domains subfamily A member 8; four-span transmembrane protein 4; membrane-spanning 4-domains, subfamily A, member 8B; membrane spanning 4-domains A8		

For research use only.

GENE DESCRIPTION	membrane-spanning 4-domains, subfamily A, member 8B
GENE ID	83661
ORIGIN SPECIES	Human
ACCESSIONS	BC022895
VECTOR	pDONR223
NCBI SYMBOL	MS4A8B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaattcgaatgactcagcaggtccggtgccaattctgtgtggtggcaccacaaatggtatcctgtgacccagggaattatgtctacgtgccctgtatccaaacagccagccgcaagtccacctagttcctgggaaccacactagtttggtgcgaatgtgaatgggcagcctgtgcagaaagctctgaaagaaggcacaaccttggggccatccagatcatcattggcctggctgcacatggcctcggctccatcatggcgacgggttctctgtaggggaatacctgtctatttattctacggaggctttccctctggggaggcttgggtttatcattcaggatctctcctgtggcagcagaaatcagccatattcttattgctgtctgtgagcagttgggttgaacatcgtcagtgcaatctgtctgcagttggagtcatactcttcacagatctaagtattccccacccatagcctaccccgactattatccttaccgtgggtgtgaacctggaatggcgatttctggcgtgctgctgtctctgctcctgagtttggcatcgcgcacatctccacttggctgccagttggtctgtgtcaatcaagcaatgtgagtgcatctatccaacatctatgcagcaaacccagtgatcacccagaaccggtgacctaccaccaagttattccagtgagatccaagcaataagtaa