

NRSN1-neurensin 1 Gene

Reference: PTXBC023514

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	p24; neurensin-1; neuro-p24; vesicular membrain protein p24; vesicular membrane protein of 24 kDa; vesicular membrane protein p24; neurensin 1		

For research use only.

GENE DESCRIPTION	neurensin 1
GENE ID	140767
ORIGIN SPECIES	Human
ACCESSIONS	BC023514
VECTOR	pDONR223
NCBI SYMBOL	NRSN1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagttcttcagcaacgtctgtgggtccaggcaggcacaggctgcagctgaggggtgtaccagcgctatggagtcgggtcctacctgcaccagtttatgaggactgtacagcctcaatttggagtaggagtagttccagatccaagatcacctaacaggtggagctcagtattct
ggaaggttggaactatctcaggtacagttttgtatcctcgggattgactgttttgccagtggtgctttctgtgcccccaaatcgaagcatttggcgaagccgattttgggtggtcgcacacatgctgtccagtttaacagtgctctggacatgtacaagctggcaggagctgtctctctgc
attggaggcagctccatggcaggggtgctgctgctgtcgtgtttgtaagagctactccaagaagaaaaattcctcagcagaagttaagaacgaatcgagacatcaaagcccacaccagccggttacaagctccagggcaggggaacaaagattccagtcactttgttc
aggggtcaaatgtccagcctactgtgcaacctga