

NOL7-nucleolar protein 7, 27kDa Gene

Reference: PTXBC023517

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C6orf90; PQBP3; RARG-1; dJ223E5.2; nucleolar protein 7; nucleolar protein 7, 27kDa; nucleolar protein of 27 kDa; polyglutamine binding protein 3; retinoic acid repressible protein		

For research use only.

GENE DESCRIPTION	nucleolar protein 7, 27kDa
GENE ID	51406
ORIGIN SPECIES	Human
ACCESSIONS	BC023517
VECTOR	pDONR223
NCBI SYMBOL	NOL7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtgcagctccgaccgcgagcgtctcgcgccccggcgctcggcggaggcgatggtgacgagggccagctggcctcggaggaggaggaggcggagcagggctgttgcctgggcagcccagcagcggcgcgccgcccggagccctggaggaagacgaggaaggggacg
atgagtttgacgatgagggcccgaggagctgactttcgccagcgcggcgggaaagcgagagaagaggagcgcgagtgccgggagaccgtgcgcagggaataaacgctcctgaaggagaagaggaaagcgcgagggagctgttcacgaacagaagaaaaaaactcctt
ccagacactatgttgagaagttaaccacagcttcacagactaacatcaagaatcgccaggaaggtgaaagaagttatgtcaaaagaaaaatgaagactgtgaaaaaggaatgactccaagaaagttaagtaaaaaagtaacagctgtcagccagaataaaagctacttgccgt
aaggctaaagaccaagatctgagagattcaaggcaacaagcagcacaagccttcatacataatcattataggccaggaaccaacagcactactgtaataagttctgtcttgcacaagaaggtaccagtgaagaagagctgctgtccagttttgaataatgcttggggaatcaa
aaaaaacaaaatgccaagaggttaaaagacggtggatggtcagaaagatgaaaactaagaagtaa