

WBP5-WW domain binding protein 5 Gene

Reference: PTXBC023544

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	WBP5; WEX6; transcription elongation factor A protein-like 9; TCEA-like protein 9; WBP-5; WW domain binding protein 1; WW domain-binding protein 5; pp21 homolog; transcription elongation factor S-II protein-like 9; transcription elongation factor A like 9		

For research use only.

GENE DESCRIPTION	WW domain binding protein 5
GENE ID	51186
ORIGIN SPECIES	Human
ACCESSIONS	BC023544
VECTOR	pDONR223
NCBI SYMBOL	WBP5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaatcctgtcaaaaatggaaggaaaaccagaaaatgagagtgaaccaaagcatgaggagagccaaagcctgaggaaaagccagaagaggagagagctagaggaggagccaaagcaaaaggaaacttttagaaaaggctgattcaatctctccaggagtttaagaaga
tatacacacaggcatttaagcaatgaagatatgttagagaagtggatgaaatagatgagataaggagagtcagaacaaacttatagtgatgcgttggagggttaatcgaaccatccttaccctatttaagttag