

FRMD3-FERM domain containing 3 Gene

Reference: PTXBC023560

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	4.1O; EPB41L4O; EPB41LO; P410; FERM domain-containing protein 3; band 4.1-like protein 4; band 4.1-like protein 4O; ovary type protein 4.1; protein 4.1O; FERM domain containing 3		

For research use only.

GENE DESCRIPTION	FERM domain containing 3
GENE ID	257019
ORIGIN SPECIES	Human
ACCESSIONS	BC023560
VECTOR	pDONR223
NCBI SYMBOL	FRMD3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcagccgaggagaccaaggggattggaataaccagatggcaaaatgcttggtaagatccagactcgacgaagcctcaattgcacatggtgaaccactgcaacagtaattgtgttctgtctgtgagactcggctccaaggtcacagctgtaacactggtgttcattgcctaaa
gaggagaacatttctgctcccttgatctccagctccccagtgaggcagccgggagatgaagatccccctagtgaggaggaagataaaataaagaagaacctttaaccatctgtaactagtgtaacccaagtgccagcctgtccccaccctgtggatgacgatgagattgacat
gctctttgactgtcctctagagcttgagttggaaagagaagacagatcatttgaggatctggaagcagatgaaacgctttttgattgctgaagaagaggagctgaaggaggctcgcctgtctttgtcgtggagctatgacattctgactggccatattcgggtgaaccactggtcaaga
gtttttccaggtccttgggtggcctgggactgctgctctttgtatttccctgctcctcctcttttggagtcagtcctcagcagtaa