

C10orf84-chromosome 10 open reading frame 84 Gene

Reference: PTXBC023577

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C10orf84; bA319I23.1; protein FAM204A; family with sequence similarity 204 member A		

For research use only.

GENE DESCRIPTION	chromosome 10 open reading frame 84
GENE ID	63877
ORIGIN SPECIES	Human
ACCESSIONS	BC023577
VECTOR	pDONR223
NCBI SYMBOL	C10orf84
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtggagtgggctgtacctctggcctaataaagtgacgtgagtcgaactcggagatgaagctacgttggaactctggacttaacttacaggaagataaaggagatgagagcatcagaaaaacagaaatcatagatttcaacagatgaacaaaaactgaaacagagtc
atgtaaatgcctatgaagagtgctctctggaattccatagatgtggaataaattcaagaattgcataaaaaacattctgaacagaaaagcacaacctcaagattcagagggaagaagaacgctccagaaaagataaattgaagaatgaaaaaatacatagtgaaacctc
aatgaaacccagtggaagagcttactcagtaatttggagtcattgacccgctgttaaaaggaaaaagttgagaagtcaggccttgaagaaggatagaccaggctgtggaggagtggaattatgagaaggctgaggaaactcagcaaccagtagctactgagagct
ggtgtaaaaattgccaaagcagttgcctgccacaacttgtaaaagccaaaaggaggtgaaattcacaggctgccgaaaaaagaagaactgcatgggggtttgaagcaagaagatgggaacaaaagcaacatgggatacatgtaa