

PSCA-prostate stem cell antigen Gene

Reference: PTXBC023582

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PRO232; prostate stem cell antigen		

For research use only.

GENE DESCRIPTION	prostate stem cell antigen
GENE ID	8000
ORIGIN SPECIES	Human
ACCESSIONS	BC023582
VECTOR	pDONR223
NCBI SYMBOL	PSCA
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaggctgtgctgttgcctgttgatggcaggtggccctgcagccagcactgccctgctgtgctactctgcaaaagccaggtgagcaacgaggactgcctgcaggtggagaactgcacccagctggggagcagtgctggacgcgcgcacccgcagtggtgcctctga
ccgtcatcagcaaaagctgcagctgaaactgcgtggatgactcacaggactactacgtggcaagaagaacatcacgtgctgtgacaccgacttgtgcaacgccagcggggccatgccctgcagccggctgccgccatccttgcgtgctccctgcactcggcctgctctctggga
cccgccagctatag