

FXN-frataxin Gene

Reference: PTXBC023633

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CyaY; FARR; FRDA; X25; frataxin, mitochondrial; Friedreich ataxia protein		

For research use only.

GENE DESCRIPTION	frataxin
GENE ID	2395
ORIGIN SPECIES	Human
ACCESSIONS	BC023633
VECTOR	pDONR223
NCBI SYMBOL	FXN
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtggactctcgggcccgcagtagccggcctctggcgtacccagccggcccagccagaccctcaccgggtccggcggcggcagagttggccccactctcgggcccgcgtggcctgcgcaccgacatcgatgcgacctgcacgccccgcccgaagttcgaac
caacgtggcctcaaccagatttggatgtcaaaaagcagagtgctctatttgatgaatttgaggaaatctggaactttggccaccaggctctctagatgagaccacctatgaagactagcagaggaaacgctggactcttagcagagttttgaagacctgcagacaagccatacagtt
tgaggactatgatgtctcttgggagtggtgtcttaactgtcaaaactgggtggagatctaggaacctatgtgatcaacaagcagacgccaacaagcaaatctggctatcttccatccagtgagactaagcgttatgactggactgggaaaaactgggtgtactccacgacggtgtcc
ctccatgagctgctggccgagagctcactaaagccttaaaaacaaactggactgtcttcttggcctattccgaaaagatgcttga