

TUSC2-tumor suppressor candidate 2 Gene

Reference: PTXBC023976

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C3orf11; FUS1; PAP; PDAP2; tumor suppressor candidate 2; PDGFA-associated protein 2; fus-1 protein; fusion 1 protein		

For research use only.

GENE DESCRIPTION	tumor suppressor candidate 2
GENE ID	11334
ORIGIN SPECIES	Human
ACCESSIONS	BC023976
VECTOR	pDONR223
NCBI SYMBOL	TUSC2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgccagcgggtccaaagctcggggcctgtggcccttcgctcggcgccggaggcgccggtcagagcagcaggagctgagcaagcttggcgccctcggcgccgagctgtgcccccttcgtattcacgcgcgcggtctatgttctatgatgaggatgggatctg
gctcacgagttctatgaggagacaatcgtcaccaagaacggcgagaagcggcgcaagctgaggcgagtgcataagaatctgattcctcaggcgatcgtgaagctggatcaccccccatccacgtggattccctgtgatcctctatgagggtga