

## TAX1BP3-Tax1 (human T-cell leukemia virus type I) binding protein 3 Gene

Reference: PTXBC023980

|                |                                                                                                                                                                                                                   |                     |                |
|----------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                                               | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                                              | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | TIP-1; TIP1; tax1-binding protein 3; Tax interaction protein 1; Tax1 (human T-cell leukemia virus type I) binding protein 3; glutaminase-interacting protein 3; tax-interacting protein 1; Tax1 binding protein 3 |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | Tax1 (human T-cell leukemia virus type I) binding protein 3                                       |
| <b>GENE ID</b>          | 30851                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC023980                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | TAX1BP3                                                                                           |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atgtcctacatcccggccagccggtcaccgctgtgcaaaagattgaaattcacaagctgcgtcaaggtgagaacttaactcctgggttcagcattggagtggaatcgaccaggaccctccagaatcccttctctgaagacaagacggacaagggtattatgtcacacgggtgtc  
tgaaggaggccctgctgaaatcgtgggtgcagattggagacaagatcatgcaggtgaacggctgggacatgaccatggtcacacagaccaggcccgcaagcggtcaccaagcgtcgaggagggtggtgctgctgctggtgacgcggcagtcgtgcagaaggccgtgcag  
cagtcctatgctgctag