

## FAM89B-family with sequence similarity 89, member B Gene

Reference: PTXBC023991

|                |                                                                                                                                                             |                     |                |
|----------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                         | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                        | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | protein FAM89B; LRAP25; MTRV; MTRV1; leucine repeat adaptor protein 25; mammary tumor virus receptor homolog 1; family with sequence similarity 89 member B |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | family with sequence similarity 89, member B                                                      |
| <b>GENE ID</b>          | 23625                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC023991                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | FAM89B                                                                                            |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atgaacgggctgccctcggcagaggcgcggcggggggggctgcgtttggccgggctccaccgctgcgcgcggcctcagcggcctcctaatgcgagcgggggctgtggcgggagctggagcgcgtctacagccagcgcagccgcgtcgtctgaacctcgaactcag  
cgctggccgcgtcgcgaaggagatgtgtctgcaggtgggggctgcggcagttggacatgtcctgtgtgcccagctgtggggcctgtacgagtcacccaggaactacaacacgtgtgccaagacgtgagctctgcccagcagctgtatcctcctccatcggacagctcctaccac  
ggatgcgggctgtctgacgagcagcgtccgatgccagcctgcctcctgaccgccacccttactgtgccccagacgacaatgccgtgaccagtgctgcaggatgccttcacatcagcctctga