

LRRC20-leucine rich repeat containing 20 Gene

Reference: PTXBC024001

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	leucine-rich repeat-containing protein 20; leucine rich repeat containing 20		

For research use only.

GENE DESCRIPTION	leucine rich repeat containing 20
GENE ID	55222
ORIGIN SPECIES	Human
ACCESSIONS	BC024001
VECTOR	pDONR223
NCBI SYMBOL	LRRC20
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

```
atgctgaagaagatgggtgagccgtggccagagtagcaaggaaagtcaacgagacgtggagagcggctctgacactctggacctggccgagtgcaagctgtctctccattggcatctacaaggtctctcggaatgtctctggccagatccacctcatcacctggctaacat
gagcttaagtcctcaccagcaagttcatgaccacattcagtcagctccgagagctccacctggagggaaacttctacaccgctccccagcgaggtcagtgccctgcagcaacctcaaggccattgacctgtcccggaaccagttccaggacttccctgagcagcttaccgacctgccc
gctgagaccatcaacctggagagagaacgagatcgtagatgtgcccgtggagaagctggccgacctgcccagcttgcgcagcatcaacctccgcttcaaccaactcaacgccgaggtgcgctgcatgccccgccgctcaagttgacatgctcatgtctccggaagcgcgaa
gagccccctacctag
```